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MANUFACTURER: XIAMEN SENYANG CO.,LTD.
Add.: No.189-1, Tongji South Road, Jimei District,
Xiamen, 361000, Fujian, P.R.China
Phone: 86-592-5028899
E-mail: susan@xmsy.net
Website: www.medicalxm.com

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NOTATKI

Anti-bedsore inflatable mattresses

Pump model: P01

Match to mattress model: S01

Do read instructions carefully before using the equipment!

Doc No. SY-RD-JS0044 Version. 1.2

Publication Date: 2023/12

INTRODUCTION

Congratulations on the purchase of your **Anti-bedsore inflatable mattresses**. The durable, high quality material used in the manufacturing of these products will ensure that you experience long-lasting and uninterrupted performance.

Intended Use

This product is intended to help and reduce the incidence of pressure ulcers while optimizing patient comfort. It also provides following purposes:

- To help and reduce the incidence of pressure ulcers while optimizing patient comfort.
- For long term home care of patients suffering from pressure ulcers.
- For pain management as prescribed by a physician.

The product can only be operated by personnel who are qualified to perform general nursing procedures and have received adequate training in knowledge of prevention and treatment of pressure ulcer.

Contraindication

Patient conditions for which the application of pressure relieving therapy on an alternation system is contraindicated are as follows:

- Cervical or skeletal traction
- Unstable spinal cord injuries

Manufacturer: XIAMEN SENYANG CO., LTD.

Address: No.189-1, Tongji South Road, Jimei District, Xiamen, 361000, Fujian, P.R. China

SUNGO Europe B.V.

Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands.

Tel: +31(0)10 3034500

Table of Contents

Symbol definition	1
Products Feature	3
Safety Precaution	4
Confirm before fixing	8
Parts name	9
Use instructions	10
Maintenance	14
Storage and disposal methods	15
Troubleshooting	16
Warranty and after-sales service	17
Specification list	19
EMC Declaration	20

Symbol definition

CE certification and notified body number

European representative information

Product production batch

Product serial number

Manufacturer information

II type equipment symbol

Mandatory reading of instructions

BF device symbol

IP21

Protected against solid foreign objects of 12.5mm and greater; protection against vertically falling water drops.

When the end user intends to discard this product, the product must be sent to the appropriate facility for recycling and recycling.

Keep away from heat and fire

Take care to prevent electric shock

Use environmentally friendly packaging materials

Environment temperature for operation or storage

Environment humidity for operation or storage

Atmospheric pressure for operation or storage

Maximum weight of the mattress load.

Do not wash the mattress.

Do not bleach the mattress.

Do not tumble dry.



Do not iron.

Do not dry clean.

1. Two-channel alternating mode

The tube mattress and the ball mattress alternate about every 6 minutes; support the body while inflating, relax the body while deflation to reduce stress.

2. Tube structure

Safety Precautions

In order to use this product safely, please read this book carefully and understand it, please be sure to follow the instructions in this manual. In order to avoid a disaster when using this product, you must know in advance the factors that may be dangerous. However, it is difficult to predict all potential dangers. Therefore, warnings about security are defined and described in this manual.

Warning If you do not follow the instructions in this symbol, you are warned that there is a risk of property damage, injury or death. Because it is especially important, it is described below as "safety precautions and warnings"

Safety warning

Warning 1 This product features a layer of soft air floating beneath the patient's body. Therefore, do not use when patients need to undergo cardiac resuscitation. When using emergency measures such as ICU, please use it under the guidance of a doctor.

Warning 2 Use this product with your doctor or specialist when using this product. In addition, if the body feels abnormal during use, or the symptoms are worse, or there is a risk of accident, stop using it immediately and consult an expert.

Warning 3 When installing the mattress, place the airway on the side of the foot; if the airway is on the side of the head, the airway may wrap around the neck and cause a serious accident

Warning 4 When using the mattress, be sure to secure the guardrail to the bed. Failure to use the guardrail may result in a fall from the bed, which may result in an accident.

Warning 5 Please do not jump on the mattress, it may cause injury or damage. Do not stand on an air mattress or knee on a mattress. There is a risk of damaging the airbag when a high pressure is applied locally.

Warning 6 When you move to a sitting position or get out of bed, be sure to ask someone else for assistance, otherwise there is a danger of falling.

Warning 7 Smoking is not allowed on the mattress and may cause a fire.

Warning 8 Please do not bundle the air outlet or press it under the mattress. When the blower tube is bent or compressed, air may not enter the air mattress and may not achieve the desired results.

Warning 9 Please do not use two or more people together. This may result in failure to perform the performance of the product and cause damage to the mattress.

Warning 10 Do not use the air pump in a damp place. Do not use water or urine to pump the pump, which may result in electric shock or malfunction. Avoid direct sunlight. In hot and humid places, keep it away from the wall or above 55cm to avoid moisture retention.

Warning 11 Please do not remove the screw of the air pump to open the casing and repair it yourself; it may cause electric shock or malfunction. In addition, the modification of the air pump without our consent may pose a serious safety hazard. Never let customers modify the air pump.

Warning 12 Be sure to plug the air pump's power plug into a European power outlet (230V 50Hz). If it is used above this voltage, it may cause a human accident or malfunction. Also, do not insert or unplug the power plug with wet hands as this may result in electric shock or a single malfunction.

Warning 13 If you do not use the air pump for a long time or need to clean the air pump, be sure to unplug the power cord from the outlet. If not it may cause accidents, electric shocks, and malfunctions.

Warning 14 Please do not pull on the power cord of the air pump, do not damage the power cord, and do not clip the power cord to the door. When removing the plug from the socket, be sure to hold the power plug and pull it out to avoid electrical accidents and fire.

Warning 15 Please put a cotton sheet on the air mattress so that the skin does not directly contact the air mattress to prevent the body fluid from contaminating the air mattress, which is not conducive to cleaning.

Confirm before fixing

Before fixing, please confirm whether the accessories are complete and whether you understand the contents and precautions of the manual.


Packing list confirmation

- ① Air pump for air mattress.....x1
- ② User Manual.....x1
- ③ Air mattress.....x1

1


2

3


Parts Name

Use instructions


Step1:

- Place the air mattress on the bed
- Connect the air tube with mattress

Step2:


- Hang the hook of the air pump on the bedside
- Connecting air pipe on pump side

Step3:

- Plug the power plug into the outlet

Warning: Dot to position device to make it difficult to operate the disconnection device

Step4:


- Turn On the power
- Adjust the Pressure

Step5:

- Adjust pressure based on patient's weight and condition
- Air pump output pressure table:

Min	1	2	3	4	5(max)
20mmHg	30mmHg	50mmHg	70mmHg	100mmHg	120- 140mmHg

Step6: Terminate operation of device

- When do not need to use the device, turn off the power switch and unplug the power plug

Warning

- Since the mattress is alternated every 6 minutes, it takes about 12 minutes for the pressure to reach the set value when a new pressure is set.
- And usually the pressure value of the mattress is at most about half of the set value of the air pump.
- Usually the higher the severity of pressure sores, the lower the pressure should be set, it is best to follow the doctor's advice, or have a professional person to operate. Direct use at maximum pressure is strictly prohibited.

Maintenance

Warning: Do not maintenance when the product is using.

Air pump output pressure check

- In the case of air pumps and air mattresses, please check the pressure once a week. In addition, if you change the location, or if you suspend operation due to power outages, etc., please do it every time, even within a week.

Mattress cleaning

- Usually use a dry rag to clean the mattress
- If necessary, wash the rag with a neutral detergent, wring the rag and clean the mattress.
- Let it dry naturally
- Maintenance is required when the patient is no longer in use.

Maintenance is required when the mattress is contaminated by the patient's body fluids.

Maintenance is required every other week when in use.

Warning:

- Dry cleaning is forbidden, otherwise it will hurt the mattress
- Please do not use a dryer to avoid deformation of the mattress.
- Do not use an iron to avoid damage to the mattress

Air pump cleaning

- Turn off the power switch and unplug the power cord.
- Usually use a dry rag to clean the air pump.

Storage and disposal methods

If you are no longer using this product, please follow the steps below to save.

1. Turn off the power switch and unplug the power cord.
2. Unplug the air tube
3. Cleaning dirty air pumps and mattresses
4. Fold the mattress and put it in the plastic bag
5. Tighten the power cord and put the air pump in a plastic bag.
6. Put the air pump and mattress in the box, together with this manual

Note:

- Please put it in a stable place, do not drop from a height.
- Please do not put heavy objects on the box.
- Please keep it in a place with low humidity.

Environmental Specifications

For transport and storage:

Ambient temperature: -20°C- +55°C

Relative humidity: 10% - 90%RH, but not requiring a water vapour partial pressure greater than 50hPa

Atmospheric pressure: 80KPa to 106KPa

For operation:

Ambient temperature: +5°C - +35°C

Relative humidity: 15% - 85%

Atmospheric pressure: 86KPa to 106KPa

Troubleshooting

Trouble	Possible reason
Do not start, do not light the power light	The power plug is not plugged into the power outlet, the power supply does not work.
Air mattress is too soft	Valve does not work Air duct off Air duct bend Adjustment knob error Damaged mattress
Air mattress is too soft	Adjustment knob error
Air pump noise is abnormally large	Put something on the air pump The air pump is in contact with other things Put the air pump on something that is easy to vibrate
Power outage	The power plug Internal fault

Warranty and after-sales service(Please read carefully)

Service life and Shelf life

- EXPECTED SERVICE LIFE OF the pump is 2 years
- EXPECTED SERVICE LIFE OF the mattress is 2 years
- Shelf life of the PVC mattress more than 10 years in the required environment.
- Shelf life of the TPU mattress is 2 years in the required environment, because TPU is an environmentally friendly material that will biodegrade after two years.

Warranty

Item	Warranty
PUMP	18 months
MATTRESS	12 months

About the repair request

First, look at the possible causes (see page 14). If you are still having trouble, turn off the product and contact your dealer.

During the warranty period:

- In the event of a malfunction in a normal use condition, your dealer, Will repair according to the provisions of the warranty.
- Please show your warranty when you are under warranty.
- In addition, if the second article of the guarantee is met, the repair is charged. Please confirm the guarantee for details.

After the warranty expires:

- Please consult with the shop or company you purchased, and if you can maintain the function through maintenance, you can repair it at a charge according to the customer's request.

Service personnel requirements

- Should have a certain understanding of the product.
- Should be able to understand the circuit schematic.
- Should be able to use screwdriver, soldering iron, and multimeter, etc.
- We (manufacturer) will provide the technical instructions, the circuit diagram, parts list to the service personnel.

Warning:

- Please don't replace the power cord by yourself and it must be replaced by professional service personnel.
- The parts should be replaced by professional service personnel when it is necessary to be replaced.
- The only safe way to turn off the power is to unplug the power cord

Specification list

Applied Part: Air mattress details description				
Mode	B01	S01	T01	T05
Inflated Size	L190xW90XH7(cm)	L186xW77xH6.5(cm)	L200xW90xH10(cm)	L200xW90xH10(cm))
Weight	1.95Kg	2.4Kg	5.31Kg	4.58Kg
Cell Material	PVC	PVC	PVC	TPU

Air Pump details description	
Mode	P01
Size	L29xW12xH10(cm)
Weight	1.15Kg
Material	ABS plastic
Power	230Vac 50Hz 12VA Power cord length: 2 meters Fuse: Time-lag type 500mA 250V Size 3.6x10mm

EMC Declaration

EMC information

Guidance and manufacturer's declaration-electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that they are used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The models device use RF energy only for their internal function. Therefore, their RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	The models P01, P05 are suitable for used in domestic establishment and in establishment directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	
<p>Warning:</p> <ol style="list-style-type: none"> 1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used. 2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. 3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. 		

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Recommended separation distances between portable and mobile RF communications equipment and the models P01, P05

The models P01, P05 are intended for use in electromagnetic environment in which radiated RF disturbances is controlled. The customer or the user of the models P01, P05 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the models P01, P05 are recommended below, according to the maximum output power of the communications equipment.

Rated maximum output Power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d=1.2 \times P^{1/2}$	80MHz to 800MHz $d=1.2 \times P^{1/2}$	800MHz to 2,5GHz $d=2.3 \times P^{1/2}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz. the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

CHANGE HISTORY

Version	Description	Date
1.0	Fist release version	2019-03